
Synicem

Cyphoplastie
Kyphoplasty

 Synimed

CE
0499

Synicem KYP

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RADIOPAQUE SURGICAL CEMENT FOR KYPHOPLASTY

Synicem KYP acrylic bone cement is recommended for the repair of vertebral bodies with structural alterations due to osteoporosis, tumour infiltration (benign or malignant) vascular malformations and selected trauma cases.

Presentation:

Synicem KYP is available in a case containing a pouch with 20 g radiopaque powder sterilized using ethylene oxide and an ampoule with 9 ml of liquid sterilized by ultrafiltration.

Advantages

- Synicem Cements are EC Marked (European Community Conformity Mark) and comply with ISO 13485 and ISO 5833 relative to acrylic cements for surgical implants.
- Synicem KYP cement combines a constant high viscosity throughout it's working time with a excellent radiopacity.
- Triple sterility protection barrier.
- Excellent tracking of cement.
- Suited for the filling of vertebral cavities generated by means of kyphoplasty balloons or similar devices.

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The Synicem KYP Système Set has been designed to perform the percutaneous treatment of vertebral bodies by means of the balloon kyphoplasty technique. The balloon kyphoplasty is used for the percutaneous treatment of vertebral body fractures (mainly in acute cases), where it is sought to restore the height of a vertebral body, or in those cases where cement application through other methods might be liable to produce a cement leakage.

We recommend the use of Synicem KYP cement for this surgical procedure.

Presentation:

The box contains the instructions for use leaflet and two Tyvek™ sealed blisters sterilized by ETO containing:

BLISTER #1

- 1 Bone puncture cannula (3 x 150 mm)
- 1 Guide wire (2 x 245 mm)
- 1 Guide wire chuck
- 1 Surgical hammer
- 1 Kyphoplasty cannula (4 x 150 mm) [working cannula]
- 1 Drill with handle (3.9 x 200 mm)
- 1 Mandrel with handle.

BLISTER #2

- 1 Kyphoplasty balloon 10mm, 15mm or 20mm
- 1 Syringe with dual scale manometer: PSI / ATM
- 1 Cement mixing bowl
- 1 Cement mixing spatula
- 3 Cement application syringes (3 ml)
- 3 Cement aspiration cannulas
- 1 Connector tube.

Advantages

- Easy operation.
- Ergonomic design.
- Easy multi-level approach: short handles do not interfere with each other.
- Full Luer thread: better attachment of components.
- Self locking mandrel.
- Exclusive guide wire optimizes needle positioning.

Ø 10mm REF 880675
Ø 15mm REF 880676
Ø 20mm REF 880677

Medical devices designer and
manufacturing specialist.

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Synergie Ingénierie Médicale S.A.R.L.